

IH 60 Nebuliser



Nebulisation of liquid medicine using compressed air technology

Rechargeable, with lithium-ion battery

Silent, small and light

For the treatment of the upper and lower respiratory tract For use in colds, asthma, respiratory disease, etc. - for fast and effective inhalation High proportion of respirable particles High nebulisation capacity (approx. 0.25 ml/min) short inhalation time Accessories: Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask, compressed air tube, mains adapter with micro-USB cable, filter, storage bag and carrying bag Can be desinfected Operating pressure: approx. 0.25-0.5 bar Medical device Li-ion battery: 2570 mAh; 3.7 V Mains operation: Input: 100-240 V~; 50-60 Hz; 0.5 A Output: 5 V === 2 A Product measurements (LxWxH): 110 x 62 x 47 mm Product weight: approx. 222 g 5 year guarantee Sales unit: 3 / Shipping carton: 4 x 3 EAN no.: 4211125 60206 2 Item no.: 602.06



german engineering IH 60

beurer