



## IH 26 Nebuliser


**german  
engineering**

**compressed  
air technology**

Nebulisation of liquid medicine using compressed air technology


**nasal  
wash**

Including nasal douche



For the treatment of the upper and lower respiratory tract

For use in colds, asthma, respiratory disease, etc.

- for fast and effective inhalation

High proportion of respirable particles

High nebulisation capacity (approx. 0.3 ml/min.)

- short inhalation time

Accessories:

Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask, compressed air tube, nasal douche, filter and storage bag

Can be disinfected

Operating pressure: approx. 0.8 to 1.45 bar

Medical device

Type: IH 25/2

Mains operation: 230V; 50 Hz; 230 VA

Product measurements (LxWxH): 166 x 141 x 148 mm

Product weight: approx. 1.4 kg

5 year guarantee

Sales unit: 2 / Shipping carton: 4 x 2

EAN no.: 4211125 60135 5

Item no.: 601.35

